


# SECU-Up Update and Options Installation Program for SECUTEST ... / SECULIFE | ST

3-348-848-29  
6/5.09



<b>1</b>	<b>General</b> .....	<b>3</b>
1.1	Software Update .....	3
1.2	Activating Options – Software Upgrade .....	3
<b>2</b>	<b>System Requirements</b> .....	<b>4</b>
2.1	SECUTEST... ..	4
<b>3</b>	<b>Installation</b> .....	<b>4</b>
<b>4</b>	<b>Connecting the SECUTEST... Test Instrument to the PC</b> .....	<b>5</b>
<b>5</b>	<b>Starting the Program</b> .....	<b>5</b>
<b>6</b>	<b>Operation</b> .....	<b>6</b>
6.1	Read Characteristics .....	6
6.1.1	Assigning an Interface Port .....	6
6.1.2	Scanning Interface Ports .....	6
6.1.3	Significance of Individual Parameters .....	6
6.2	Run Update, Change Language .....	7
6.3	Installing Options .....	7
6.4	Report Options with the Optional DA-II Printer Adapter .....	7
6.5	Configuration .....	8
6.5.1	Setting the Clock .....	8
6.5.2	Assigning Names .....	8
6.5.3	Limit Values .....	8
6.5.4	Miscellaneous .....	8
6.6	Database (DBmed option) (only useful for <b>SECUTEST SIII (+) / SECULIFE ST</b> ) .....	9
6.6.1	Uploading Test Specifications .....	9
6.6.2	Reading Out Test Specifications .....	10
6.6.3	Deleting the Database .....	10
6.7	Detection of probe at protective conductor (SK5 option) .....	10
6.8	Direct Print-Out (SECU-dd option) .....	10
6.9	Resetting the SECUTEST... ..	10
6.10	Exit the Program .....	10
6.11	Help Function – Online Manual .....	10
<b>7</b>	<b>Product Support</b> .....	<b>11</b>

# 1 General

The modern design of the SECUTEST.../ SECULIFE ST test instrument allows for updating of device software to the latest standards and test specifications. Customer suggestions result in continuous improvements to test instrument software as well.

## 1.1 Software Update

This program allows for quick, on-site updating of your test instrument software. Currently activated options are updated as well. You can choose to perform a complete update for all software in your SECUTEST... test instrument, or update only individual modules: sequences, schematic diagrams, on-line help, balancing, operating system.



### Note

We recommend the automatic update sequence which updates all modules, one after the other. If, for example, the balancing module is selected for updating, an update must also be run for the operating system module. If you forget to run this second update, the following error message appears: “command error”.

## 1.2 Activating Options – Software Upgrade

The options listed below can be activated after the appropriate password has been entered. The desired options must be ordered from GMC-I Messtechnik GmbH with the enclosed registration, upon which we forward the required passwords to you.

The software allows for activation of the following options:

- **Test Sequence IEC 601 (IEC 601 option – article no. Z853G)**  
Enables measurements in accordance with IEC 601 with the SECUTEST®0751/601S and as of SECUTEST®SIII test instruments. The test sequence in accordance with IEC 601 is activated.
- **Database (DBmed option – article no. Z853H)**  
Test sequences can be configured on-site and performed in accordance with the respective requirements within a given selector switch position. Configurations for various test sequences are stored to the test instrument and can be activated at any time in the future. Measurement values acquired during the test sequence are also stored to the test instrument. These values can then be printed out onto report form templates which are stored at the test instrument.
- **Detection of probe at protective conductor (remote control) (SK5 option – article no. Z745K)**  
Protective conductor measurement is expanded to include the function: “automatic measuring point recognition”. This function is quite useful if several protective conductor connections need to be tested without having to manually restart the measurement for each measuring point. The test instrument recognizes whether or not the probe is in contact with the protective conductor during protective conductor measurement, and indicates the two possible conditions with different acoustic signals. The test sequence runs automatically.  
Rapid repetitive signal: probe in contact with protective conductor.  
Slow repetitive signal: measurement completed, next measuring point.
- **Direct Print-Out (SECU-dd option – article no. Z853L)**  
Test results are read out directly via the RS 232 interface after each test (after individual test, or at the end of a test sequence).

## 2 System Requirements

### Software

- MS WINDOWS 2000 or XP.

### Hardware

- IBM compatible Windows PC, Pentium processor with 200 MHz or better and at least 64 MB RAM
- SVGA monitor
- Hard disc with a least 20 MB available memory
- Microsoft compatible mouse
- Connection to internet or CD-ROM drive

### 2.1 SECUTEST...

- No restrictions
- Not for SECUTEST S2N+ und S2N+10
- Not for SECUTEST -DC, -med, 0701 und ... 0701/0702S

## 3 Installation

---



### Note

The latest version of the device software is available on our homepage.

---

- ⇨ Extract the zip file
- ⇨ Close all running applications before starting the installation .
- ⇨ Run setup.exe.
- ⇨ Follow the installation guide.
- ⇨ Open the **README** file which contains the latest information which may not have been available when these instructions were printed. Use the MS WINDOWS Notepad.

A setup window appears with a recommended directory to which the program is automatically installed, if you click on the **Continue** button.

You can also enter a different directory or disc drive here, and then start automatic installation.

A message appears after the program has been successfully installed, which must be acknowledged with **OK**.

After installation is complete, a program group is automatically added to the “WINDOWS Program Manager” which is identified with an icon which depicts a key-operated switch.

## 4 Connecting the SECUTEST... Test Instrument to the PC

After the program has been successfully installed, connect the test instrument to the PC as follows:

- ⇒ Switch the PC and the test instrument off.
- ⇒ Connect the PC to the SECUTEST... with the Z3241 serial interface cable.
- ⇒ Secure both interface plugs.
- ⇒ Switch the PC and the test instrument on.

## 5 Starting the Program

WINDOWS offers three different possibilities for trouble-free starting of the program:

- ⇒ Start the program from the group window which was added to the program manager during installation by double clicking on the **program icon**.
- ⇒ Start the program from the corresponding directory in the file manager by double clicking on the **Z853G.EXE** program file.
- ⇒ Start the program with the **RUN** command in the "FILE" menu in the program manager or the file manager.

See your MS WINDOWS manual for detailed instruction concerning the three ways in which applications can be started.

After the program has been started, a window with the following message appears: "Please connect the test instrument to a free COM port".

If the PC and the test instrument have already been connected and connection has been acknowledged, the initial application window shown in chapter 6.1 appears. There are no entries in the information fields at first, and the options menus do not appear.

First select a user interface language for the update program in the "Program" menu under **Language**, *Deutsch, English, Française, Italiano, ...*

## 6 Operation

### 6.1 Read Characteristics

Settings

Program Update Install Options Test results configure Help

Tester

Serial Number

Name

Options

COM 1

Read settings

Variant

customer

revision

Sequence

Help

Picture

Abgleich

System

#### 6.1.1 Assigning an Interface Port

⇒ If you do not know which port your test instrument is connected to, proceed as described in chapter 6.1.2.

Select the port to which your test instrument is connected, i.e. COM 1, 2, 3 or 4. The default setting is COM 2. Activate the **read characteristics** button. If the test instrument has been successfully connected to the PC, entries now appear in the information fields.

#### 6.1.2 Scanning Interface Ports

⇒ Activate the **read characteristics** button.

The program scans all ports to which a test instrument of the type SECUTEST... might be connected. Each unoccupied port generates an error message which must be acknowledged. If connection to the test instrument is successfully established, entries appear in the information fields.

#### 6.1.3 Significance of Individual Parameters

Serial Number	Test instrument serial number
Name	See chapter 6.5.2 on page 8
Options	Previously installed and activated options
Sequence	Date of last installed test sequence
Help	Date of last installed on-line help
Diagrams	Date of last installed schematic diagrams
System	Date of last installed operating system
Balance	Date on which last balancing was performed
Variants	Test instrument variants, e.g. instrument with 25 A test current

## 6.2 Run Update, Change Language

---



### Attention!

Neither the test instrument nor the PC may be disconnected from the electrical power supply during uploading (“Update” / “Install options”).  
All other WINDOWS programs must be closed prior to updating.

---

Select either **automatic** in the “Update” menu, or a specific module which you would like to activate.

If automatic is selected, automatic updating of all modules, one after the other, is performed. If no updating is required for individual modules, or if updating has already been completed, the following message appears:

*“Module ‘XYZ’ is current”.*

If you would like to install a different language for the test instrument user interface, the desired language can be selected in the “Update” menu under **SECUTEST language**, *Deutsch, English, Français, Italiano, ...* and uploaded to the test instrument.

## 6.3 Installing Options

An automatic update should be run before options are installed (activated), see chapter 6.2. After the update has been completed, select the desired option from the “Install options” menu.

An entry field appears into which the appropriate password must be entered.

The following message appears:

*“Activating ‘XXX’ option”.*

## 6.4 Report Options with the Optional DA-II Printer Adapter



### Using the DA-II

Test instruments which are not equipped with a parallel port can be connected to commercially available Centronics printers with the help of the DA-II printer adapter (Z745M), allowing for on-site print-out of test reports and more.

This is only possible with test instruments which have been designed for direct report generation. All SECUTEST... test instruments as of manufacturing revision level 8/99 fulfill this prerequisite.

## Uploading, Editing and Saving Report Templates

The „Test result“ menu allows for storage of up to 4 **report templates** in the SECUTEST... Report templates are uploaded from a file at a PC to the SECUTEST... (“Load file”). Templates can also be uploaded from the SECUTEST... to a PC (“Templates from Secutest”), and then edited and re-saved.

After completion of a test (“Passed” or “Failed” appears at the instrument’s LCD), the test results are read out to one of these report templates in the form of **report data** (depending upon the selector switch position or test regulation) via the RS 232 interface. The report menu in the SECUTEST... is activated with the  and  keys to this end. Report data can only be read out to a Centronics printer or a PC if a DA-II adapter has been connected to the RS 232 interface. The report can, for example, be displayed with the help of a terminal program such as Hyperterm for WIN NT.

## Reading Out, Printing and Saving Test Results / Test Data from the (P)SI-Module

Test results saved to the SECUTEST... can be displayed, edited, printed or saved after reading them into a PC using the desired report template.

Data can be processed immediately after completion of the test, or from the database (optional database in the SECUTEST...). Data can be saved with or without a report template (e.g. for further processing with PS3).

Test results saved to the (P)SI-Module can also be read out, printed, saved or processed by means of a report template.

Easy method of printing out reports:

Activate the „test result“ menu in the update and options installation program. When the test result is displayed in the SECUTEST... select the print function (press key „Cursor up“, then set the cursor on „print“ and press ENTER).

## 6.5 Configuration

### 6.5.1 Setting the Clock

Date and time can be set here in the following format:

dd.mm.yy; hh:mm (d = day, m = month, y = year; h = hour, m = minute).

### 6.5.2 Assigning Names

In the “Assign names” menu you can, for example, enter the name of the inspector who normally works with this test instrument, or the name of the department which is responsible for testing.

The name may not exceed a length of 8 alphanumeric characters including spaces.

The name appears as the first entry in reports which are transmitted to the PC.

### 6.5.3 Limit Values

The limit values for the internal test items „short-circuit test“ and „voltage reference“ can be modified here.

### 6.5.4 Miscellaneous

„Use designation“

When starting a test, the user is requested to enter a designation. This designation will then be saved and printed out in the test report (in addition to the ID no.)



## 6.6 Database (DBmed option) (only useful for SECUTEST SIII (+) / SECULIFE | ST)

If the **database** option has been activated, the “Database” menu can be opened. The limit values stored to the test instrument are read out when the “Database” menu is selected. Select the appropriate test specification, safety class and application part type, in order to display the desired limit values. Measurements for which no limit values have been stored to memory will not be performed.

The screenshot shows the 'Database' software interface. The title bar reads 'Database' and the menu bar includes 'Exit', 'Store sequence', 'Read sequences', and 'Clear data'. The main window is titled 'classify appliance' and displays 'IEC 601' in a dropdown menu. Below this, there is a 'general' section with a 'Class' selector set to 'I' and a '0 Applied part's' section with radio buttons for '0 II', '0 B', '0 BF', '0 CF', and a checkbox for 'R-ISO'. On the right, a 'Sequences stored' table lists five entries with their respective measured times and dates.

Sequences stored	
1: measured:	22:17:56 11.12.00
2: measured:	05:44:11 12.12.00
3: measured:	05:44:29 12.12.00
4: measured:	21:10:57 12.12.00
5: measured:	03:00:25 13.12.00

Parameter	Value		
Protective Earth Resistance	<n. 300Ω		
Protective Earth Current	<0.500mA		
Insulation Resistance	>2.000MΩ	Protective Earth Current SFC	<1.000mA
Insulation Voltage	>0500V	Enclosure Leakage Current	<0.100mA
Residual current	<3.500mA	Enclosure Leakage Current SFC	<0.500mA
U-AC		Patient Leakage Current AC	<0.100mA
U-DC		Patient Leakage Current AC SFC	<0.500mA
Equivalent Leakage Current		Patient Leakage Current DC	<0.010mA
mains on applied part	<05.00mA	Patient Leakage Current DC SFC	<0.050mA
High voltage		Patient Auxiliary Current AC	<0.100mA
		Patient Auxiliary Current AC SFC	<0.500mA
		Patient Auxiliary Current DC	<0.010mA
		Patient Auxiliary Current DC SFC	<0.050mA

Test results (measured) are also included in the list of test regulations and test specifications. The test results are contained in report templates which can be edited via the „Test result“ menu. Branching is established to the „Test result“ menu by double clicking a report template.

### 6.6.1 Uploading Test Specifications

Create an individualized test sequence by entering or changing the following parameters in the indicated order:

1. Test Specification, 2. Safety Class, 3. Application Part Type, 4. Limit Values

When the **Upload test specification** command is selected you are requested to enter an ID number. This allows for unambiguous identification of the test sequence. The ID number can be entered as a series of characters containing numbers and/or letters. The default setting uses consecutive numbers with which the data are stored to the test instrument.



#### Note

The selector switch position must be in agreement with the test specification when the stored test sequence is run.

## 6.6.2 Reading Out Test Specifications

The limit values stored to the test instrument are displayed when the respective ID number is selected in the “Stored test sequences” field with the cursor.

---



### Note

If the  $\Omega$  symbol is not displayed, the appropriate font has not been installed. Install the **linedraw** font to the control panel.

---

## 6.6.3 Deleting the Database

---



### Attention!

If the **Delete database** button is activated, all test sequences which have been stored to the test instrument are deleted, i.e. limit values as well as measured data. Limit values in accordance with DIN VDE remain.

---

## 6.7 Detection of probe at protective conductor (SK5 option)

Upload the remote control module to the test instrument as follows:

- ⇒ Select **Remote control** in the “Options” menu.

Expand protective conductor measurement to include the function “automatic measuring point recognition” as follows:

- ⇒ Activate the “Auto measuring point” operating parameter in the corresponding configuration page.

## 6.8 Direct Print-Out (SECU-dd option)

Test results are read out directly via the RS 232 interface after each test (after individual test, or at the end of a test sequence). In contrast to the result of a test sequence where only the lowest value of each test is indicated, every measured value is documented in this case.

## 6.9 Resetting the SECUTEST...

If unexpected problems should occur during updating or options activation, the operating system can be restarted. This results in the same status which occurs after operating voltage has been applied to the test instrument. Select **Reset** in the “Program” menu.

## 6.10 Exit the Program

Select **Exit** in the “Program” menu.

## 6.11 Help Function – Online Manual

If a printed manual is not available on site, the online manual can be queried.

- ⇒ Select item **Contents** in the „Help“ menu.

## **7 Product Support**

When you need support, please contact:

GMC-I Messtechnik GmbH  
Product Support Hotline  
Phone +49 911 8602-112  
Fax +49 911 8602-709  
E-Mail [support@gossenmetrawatt.com](mailto:support@gossenmetrawatt.com)

---

Printed in Germany • Subject to change without notice

 **GOSSEN METRAWATT**  
GMC-I Messtechnik GmbH  
Südwestpark 15  
90449 Nürnberg • Germany

Phone +49 911 8602-111  
Fax +49 911 8602-777  
E-Mail [info@gossenmetrawatt.com](mailto:info@gossenmetrawatt.com)  
[www.gossenmetrawatt.com](http://www.gossenmetrawatt.com)