

Brief instruction for SQL database installation

1. System requirements:

Recommended: Windows 7 or higher with .NET Framework 4.0

2. Installation SQL Server R2 Express:


Subsequently, the download link for the SQL Server 2008 R2 Express. Depending on the operating system (x86 or x64) a respective selection is required: SQLEXPRTW_x86.exe or SQLEXPRTW_x64.exe.

<http://www.microsoft.com/en-us/download/details.aspx?id=30438>

Complete installation instructions:

[https://technet.microsoft.com/en-us/library/ms143219\(v=sql.105\).aspx](https://technet.microsoft.com/en-us/library/ms143219(v=sql.105).aspx)

- a) Select following installer file: ENU\x64 \SQLEXPRTW_x64_ENU.exe“ (WT= With Tools).




Next steps after successful download:

- Double click .exe-File
- Start of installation process


- c) Installation window for SQL Express: Double-click on “New SQL Serber stand-alone installation or add features to an existing installation”.




- a) Start and execution of the actual installation




d) Acceptance of license terms




e) Feature Selection. Select All!







f) Selection of „Named instance“, Enter: „SQLExpress“




g) Leave all default settings as they are in the next steps.







a) Completion of installation process



Instruction for installation of the SMARTCOLLECT software please see „SMARTCOLLECT Manual“