

DB-COMFORT UPGRADE

FEATURE PACK FOR
SECUTEST AND SECLIFE

(FROM FIRMWARE 2.1)



SOFTWARE UPGRADE FOR SECUTEST AND SECLIFE TEST INSTRUMENTS

Main Features of Our Test Instruments

	SECUTEST			SECULIFE		
	BASE	PRO	PRO BT COMFORT	ST BASE	ST BASE 25	ST PRO
	M705B	M705C	M705E	M693F	M693G	M694C
10 A R _{PE} test current	—	■	■	■	—	■
25 A R _{PE} test current	—	—	—	—	■	—
Touchscreen keyboard	—	■	■	■	■	■
Second test probe	—	■	■	■	■	■
Voltage measuring inputs	—	■	■	■	■	■
Database expansion	—	■	■	■	■	■
Antimicrobial properties	—	—	—	■	■	■
DB-Comfort	□	□	■	□	■	■
Bluetooth	□	□	■	□	□	□

□ optionally available

Special Features and Performance Characteristics of DB-Comfort Feature Pack

■ NEW ADDITIONAL “MEDICAL” DATABASE OBJECT

Device with extended entry options

■ USER-DEFINED TEST SEQUENCES

The number of user-defined sequences has been increased to 24.

■ MOVING TEST OBJECTS

Moving (medical) devices within the tree can be initiated by pressing and holding the respective element in the tree display.

■ TOUCH-EDIT

“Editing” of a (medical) device can be started by pressing and holding the respective element of the detail display in the main window.

■ AUTO-STORE

The auto-store function can be activated in the setup menu so that the results of automatic testing are saved immediately under the selected test object.

■ PUSH PRINT

A computer connected to the test instrument can cause the SECUTEST to enter a mode in which data are transmitted directly to the connected PC instead of being stored at the tester.

■ QUICK EDIT

The quick edit option can be activated when entering a new test object so that all other fields can be filled out immediately after entering the ID number.

■ SEARCH

Searches started with the “Search All” softkey now search the new “UDI” field (unique device identification) at medical devices as well.

GMC-I Messtechnik GmbH

Südwestpark 15 ■ 90449 Nürnberg ■ Germany

Phone +49 911 8602-999 ■ Fax +49 911 8602-125